



1 Scope

Non-sterile, reusable product (up to 5x)
Monopolar HF Arthroscopic Electrodes, Article No.: 8M910331 (833-1), 8M910332 (832-1), 8M910333 (833-4), 8M910334 (833-6), 8M910335 (832-2), 8M910821 (802-8), 8M910822 (802-6), 8M910833 (805-6), 8M910834 (805-1), 8M910903 (802-9), 8M910913 (802-1), 8M910914 (800-4), 8M910915 (800-5), 8M910916 (802-5), 8M910917 (802-7), 8M910918 (802-7), 8M910919 (802-7), 8M910924 (802-7), 8M910945 (805-1), 8M910954 (800-2), 8M910971 (833-2), 8M910976 (805-5), 8M910978 (800-1), 8M911039 (805-2).
Bipolar HF Arthroscopic Electrodes, Article No.: 8M910808 (879-2), 8M910809 (879-3), 8M910843 (872-4), 8M910904 (870-3), 8M910905 (870-9), 8M910906 (879-1), 8M910983 (879-1), 8M910984 (879-9), 8M910985 (879-1), 8M911008 (870-1).
Maximum rated voltage of accessory (Umax):

Article No.	Umax/Pmax	Note
8M910331 (833-1), 8M910332 (832-1), 8M910333 (833-4), 8M910334 (833-6), 8M910821 (802-8), 8M910822 (802-6), 8M910833 (805-6), 8M910834 (805-1), 8M910903 (802-9), 8M910913 (802-1), 8M910914 (800-4), 8M910915 (800-5), 8M910916 (802-5), 8M910917 (802-7), 8M910918 (802-7), 8M910919 (802-7), 8M910924 (802-7), 8M910945 (805-1), 8M910954 (800-2), 8M910971 (833-2), 8M910976 (805-5), 8M910978 (800-1), 8M911039 (805-2)	1,3 kVp/ 120 W	
8M910808 (879-2), 8M910809 (879-3), 8M910843 (872-4), 8M910904 (870-3), 8M910905 (870-9), 8M910906 (879-1), 8M910983 (879-1), 8M910984 (879-9), 8M910985 (879-1), 8M911008 (870-1)	4,3 kVp/ 120 W	
8M910904 (870-3), 8M910905 (870-9), 8M911008 (870-1)	1,0 kVp/ 120 W	CAUTION Not to be used for vaporization, only for cutting.
8M910808 (879-2), 8M910809 (879-3), 8M910843 (872-4), 8M910904 (870-3), 8M910905 (870-9), 8M910906 (879-1), 8M910983 (879-1), 8M910984 (879-9), 8M910985 (879-1), 8M911008 (870-1)	1,0 kVp/ 180 W	

See also label or catalogue
In any combination with another electro-surgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used. (See also section "Intended Use")
See these instructions for use, the label, or the current product catalog for the maximum rated voltage of the product.

If there is any uncertainty, contact the manufacturer.
Before use, read the entire IFU of this product and any user accessory as well as HF generator.
All requirements, safety notices, and warnings included in the respective IFUs have to be followed strictly.
This medium is not intended for users in the USA.

2 Intended Use

Product intended for application in arthroscopic surgery.
For cutting, coagulating and vaporizing (exceptions see table above) of biological tissue, while applying irrigation fluid (see below).

Only to be used by skilled medical professionals who have been introduced to the product, e.g., using this IFU.
Connect the product to the monopolar/bipolar output of the HF generator via an HF cable, HF handle, or handle with suction.
Activation takes place using a foot switch or an HF handle.
It is strictly prohibited to manipulate the electrode tip.

Indication:
Product for arthroscopic application for cutting, coagulating and vaporizing biological tissue.
The product is intended to be used in large joints, e.g., knee, hip, shoulder, etc.

Contraindication:
The product is not intended for direct contact with the heart, central circulatory system or central nervous system and the spine area.

Use the product without the application of conducting irrigation fluid is strictly prohibited.
Regarding products with suction, the suction system always has to be connected and activated. The product is not intended for use in small joints.

The application of high-frequency current may interfere with cardiac pacemakers and in vivo heart defibrillators, so affected patients must consult a cardiologist before the intervention.
If arthroscopic intervention is contraindicated for the patient, he/she is not allowed to use the product. The patient must NOT come into contact with grounded metal objects during usage.

Side effects:
Electrosurgical interventions may lead to undesired damages or undesired burns of surrounding tissue.
Tissue damage due to extensive forces during blunt dissection is possible.

It is the sole responsibility of the physician to include this aspect in their professional decision regarding the treatment process and the patient's consent.

All electro-surgical products potentially may lead to muscle stimulation during the application.
The product has been designed in such a way that the risk of this undesirable effect is minimized, however, muscle stimulation still can lead to an unexpected movement of the patient in the surgical field.

It cannot be excluded, that the electrode tip wears off earlier than expected when contact with bony surfaces is done intensively.

Connection and activation:
Before using monopolar products, ensure that a designated HF neutral electrode is correctly applied via the instructions, safety and warning notes in the instructions for use of the accessories and the HF generator used must be followed.
The product is connected to the corresponding monopolar/bipolar output of the HF generator via a suitable monopolar/bipolar HF cable or a suitable monopolar/bipolar HF handle.
Activation is via the buttons on the HF handle or via a foot switch on the HF generator.

Combination / Compatibility:
Before use, the compatibility of the product with the intended HF handle/HF cable and HF generator has to be verified.
In case of uncertainties, contact the manufacturer of this product or the manufacturer of the used accessory or HF generator.

The frequency of the used HF generator shall not exceed the maximum frequency of 4 MHz and not exceed the maximum rated voltage of the accessory (see section "Scope").
The product can be connected to the following combination products of the manufacturer REGER, or to a comparable combination product: Article No. 800-0 to 800-10, 801-0 to 801-10, 832-1, 832-5, 833-1 to 833-9 are only to be used with REGER HF cable Article No. 92080 to 92083, 92181.

Feasible HF generators for the aforementioned monopolar articles:

ERBE	KLS Martin	EMED	BOWA	EMED ARTo
Covidien	ValleyLab	Tekno	Berchold	EMED ARTo distributed by REGER

Feasible HF-generators for aforementioned bipolar articles:

ERBE VIO	KLS Martin	EMED	BOWA	EMED ARTo
Depuy Mitek VAPR	Arthrocare	Tekno	Berchold	EMED ARTo distributed by REGER

Only use in combination with compatible products and connectors.
The instructions, safety and warning notes in the instructions for use of the accessories and the HF generator used must be followed.
Follow instructions and warnings in the IFU of HF cable and HF generator.

IMPORTANT:
It is recommended to use a smoke evacuation system.

Hand with utmost care.
This does not only apply for the duration of the surgery but also for the complete duration of storage, processing, and transport as well as during the process of connecting the product with the respective HF handle, HF cable, and/or HF generator.

This applies especially to the thin tips and other sensitive areas, e.g., the insulation.
Improper use will result in immediate loss of warranty.
No liability is accepted for any damage that may occur.

3 Safety Notices – WARNING!

Make sure the User Manual of the HF generator has been read carefully before using the product.
Proceed according to the maximum rated voltage stated in this IFU or label.
Before each use, the product must be cleaned, disinfected, and sterilized according to a validated procedure (DIN EN ISO 17665). (See section "Reprocessing: Cleaning, disinfection, and sterilization".)

In case of uncertainties, contact the manufacturer.
Before usage, a visual inspection and function test has to be done (see section "Visual Inspection and Function Test").

Before usage, a visual inspection has to be done for pressure points of damage on the product.
In case of damage, deformation, or similar is detected on the product and its insulation, the product must not be used.
It has to be replaced with a new one.

Irrigation fluid for monopolar products:
During use, the active electrode must always be completely surrounded by 100% electrolyte-free rinsing solution, e.g., use Purisolve. Do not use without rinsing solution.

Irrigation fluid for bipolar products:
When in use, both the active electrode and the neutral electrode of the product must always be completely surrounded by 100% conductive flushing liquid.
Use 0.9% saline solution as flushing liquid, for example.
NEVER use non-conductive irrigation fluid (e.g., sterile water, Glycine, Purisolve, etc.) for bipolar products.

Ensure sufficient and adequate flow of irrigation fluid to the joint (at least 100 ml/min).
At least one (1) sterile backup product has to be available for the application.

It is the responsibility of the user to determine the appropriate product size according to their professional judgment and based on the patient's specific indication, preferred surgical technique, history, etc.

Only after application, verify that the product is inserted firmly in the HF handle / HF cable.
This must be done carefully, to avoid damages to the product and/or injuries to the patient, surgical personnel, or third party.

Excessive force can damage the product.
Therefore, the product has to be observed during the complete application.
Exclusion:
Do not activate the product as long as it is in contact with metal objects and/or optics.
During an electro-surgical intervention, the patient must not come into contact with grounded metal objects such as surgical desk frames, instrument trays, etc.

Make sure that there are no flammable or combustible substances (e.g., disinfectants, gases, etc.) in the immediate vicinity during the entire application, as heating or sparking could cause a fire or explosion.
Do not touch the product tip as long as it is activated. After switching off the electro-surgical current, the product tip still may be hot, so this could lead to undesirable burns. While the product is activated, do not insert or withdraw the tip from the surgical site. Inadvertent activation of the product or moving the tip outside of the field of vision can lead to unintentional injury of tissue.

During activation, make sure the product does not contact the arthroscopic as thereby the product quality or arthroscopic solution may be damaged.
Do not activate the product uninterrupted over a longer period.
Do not use metal cannulas as they may damage the product insulation or create an alternative current path (capacitive coupling), resulting in unintentional burns.

It is recommended to use plastic cannula systems.
It is the responsibility of the user to select the lowest possible power setting on the HF generator to achieve the desired effect during the respective operation.
During the surgical procedure, mechanical forces may cause deformation or wear of the product.
The use of the product is disallowed during the application, the product has to be exchanged for a new one.
Abrasion, adhesion of tissue, discoloration, sooting, etc. do not represent a reason for complaint and do not permit claiming the manufacturer's warranty.

4.1. Maximum number of reprocessing cycles
Due to the materials used and due to the intended use of the products, the articles listed in the Scope of this IFU only are allowed to run through a maximum of up to 5 reprocessing cycles (refer to section "Visual Inspection and Function Test").
When applied according to the Intended Use, the product undergoes natural wear and tear, considering manner and duration of the application as well as manner and frequency of reprocessing.
Therefore, a visual inspection and function test has to be done before each use. (Refer to section "Visual Inspection and Function Test")
The visual inspection and function test, in particular the condition of the insulation and the product tip, is decisive for whether the product may be used.

4.2. Time requirements for cleaning and disinfection.
Preparation for cleaning, pre-cleaning and automated cleaning and disinfection has to be done immediately after the application, however not later than 1 hour after the application. Avoid idle time.

4.3. After the Application
Clean and disinfect product immediately after the application. However not later than 1 h after application. After application of the product, deposit the product carefully (protection of lifetime of the product). After application, separate contaminated product and deposit it in a suitable container (deposit means "do not drop"). Immediately remove gross stain. Immediately mark damaged or defective products. Accessory, that do not fit for the sieves of the cleaning and disinfection device (CDD), shall be deposit separately in suitable containers. Close firmly all disposal-transport containers, in order to avoid drying of stain. Organize transport of contaminated products that contact contamination of the transport ways and environment also is avoided (closed transport). Unused reusable products have to be cleaned as well. Refer to ARK (German organization for instrument reprocessing). Red Brochure, pages 30-32. Take care, all and any transport containers are cleaned and disinfected after transport as well.

4.4. Validation of (Re)Processing
The following validated processing procedure is recommended.
Equivalent deviant processes are possible.
Then, it is the sole responsibility of the user to safeguard the suitability of the applied procedure by suitable means (e.g., validation, routine monitoring, verification of material compatibility, etc.).
Automated disinfection always is preferable to manual disinfection.
The following procedure was validated according to EN ISO 17665.

Additional applicable processing requirements specific to the respective clinical place (operator) as well as national or country-specific regulations have to be followed as well.
Never use sharp objects for cleaning.
Disinfectants always have to be rinsed and removed carefully.

4.5. Preparation for Cleaning
Remove the product from its packaging.
Place it in a container provided for cleaning.
It is not necessary to disassemble the product.

4.6. Pre-Cleaning
Pre-clean the product immediately. At the latest, however, 1 h after completion of the application.
Use tap water (potable water quality) (<40°C) and aldehyde-free, non-fixing disinfectants if applicable.
Thoroughly remove the surface stain with a soft brush or synthetic-free, as otherwise particles or dried secretions may adhere.

Ensure that hard-to-reach areas are well cleaned and rinsed several times.
This may make subsequent cleaning and sterilization difficult or impossible.
Drain Cavities and lumen have to be rinsed intensively using at least 3x20 ml cold tap water (<40°C) with the aid of a rinsing adapter (e.g., from the company Medisafe), or with a syringe or with a water jet pistol (>30 Sec.).

This pre-cleaning step always has to be done prior to the manual cleaning or cleaning with the cleaning and disinfection device (CDD).

4.7. Manual Cleaning and Disinfection
Prepare an immersion bath with a suitable fluid cleaning agent.
Use a cleaning agent that is compatible with the disinfectant and suitable for immersion baths.
Only use agents suitable for medical devices made from metal and plastics with a pH value between 5.5 and 12.3.

Recommendation: Cleaning agent gazygarm® (Schülke & Mayr) and disinfectant Korsolve Plus.
Do not use high alkaline cleaning agents.
These will impair the lifetime of the product.

• Prepare immersion bath with cleaning agent according to the specific cleaning agent IFU.
• Prepare a separate immersion bath with disinfectant according to the specific disinfectant IFU.
• Immerse the product completely in the ultrasonic bath with a cleaning agent (e.g., 0.5% gazygarm®).

• Clean product in an ultrasonic bath using a sonication time of 5 Min. and a frequency of 45 kHz.
• Follow all instructions outlined in the IFU of the cleaning agent, disinfectant, and ultrasonic bath.
• Ensure the product will not touch other products or parts in the ultrasonic bath.
• Ensure sonic shadows in the ultrasonic bath are avoided.
• Then, clean the product with a soft brush under cold running tap water (<40°C).

• Intensively rinse cavities and lumen with a water jet pistol (>30 Sec.) or similar for at least 1 Min.
• Afterwards, rinse the product thoroughly for at least 1 min. under running tap water (>40°C) to remove any residues of the cleaning agent.
• Inspect the product visually on the remaining stain.
• In case stain is still present, repeat the aforementioned cleaning steps as long as it needs until no soil is present.

• Afterwards: immerse product completely in a disinfectant bath including e.g., Korsolve Plus 3%, for at least 15 Min.
• Follow the manufacturer's data for residence time.
• Ensure the disinfectant will contact all areas on the product.
• Rinse cavities and lumen several times, which means at least 3x with 20 ml each of disinfectant bath fluid.

• Afterwards: Rinse the product thoroughly for at least 1 Min. with demineralized cold water, to remove all disinfectant residues.
• Additionally: Rinse all narrow and areas difficult to access, all cavities and lumen with a syringe several times (at least 3x) using each time 20 ml cold demineralized water.

• Dry product with a lint-free wipe and sterile compressed air.
• Dry cavities, lumen, and channels with sterile compressed air.

4.8. Automated Cleaning and Disinfection:
Only use cleaning and disinfection devices (CDD) with proven efficiency according to EN ISO 15883. Follow the data of the manufacturer of the cleaning and disinfection machine.
Only use agents suitable for medical devices made from metal and plastics with a pH value between 5.5 and 12.3.

Recommendation: neodisher® mediclean forte (Dr. Weigert GmbH & Co. KG).
Apply program for thermal disinfection.
Follow instructions and data regarding program choices and machines.

Product has to be stored safely and protected against mechanical damages during automated cleaning and disinfection.
Do not clean together with sharp-edged or pointed objects.
Deposit product in a suitable rinsing basket.

Follow data for loading of the cleaning and disinfection device (CDD).
Use rinsing adapters for products with the lumen and connect them according to the instructions in the User Manual of the cleaning and disinfection device (CDD).

Cleaning Program
Start the program course with the following parameters:

- 1 Min. pre-rinsing with cold water
- 3 Min. pre-rinsing with cold water
- Emptying
- 5 Min. cleaning at 55°C with 0.5% alkaline cleaning agent
- Emptying
- 3 Min. neutralization with warm tap water (>40°C) and neutralizer (0.1% Neodisher® Z)
- Emptying
- 2 Min. interim rinsing with warm demineralized water (>40°C)

Disinfection Program
Automated thermal disinfection considering national requirements regarding AO value (see EN ISO 15883, AO value >3000).
5 Min. cleaning at 92°C +/-2°C

Drying
• 30 Min. at 90°C
• Remove rinsing adapter
At the end of the program course, remove the product and inspect it on the remaining stain.
Immediately after removal of the product and immediately after additional drying in a clean place, put the product in a single-use sterilization packaging (double packaging) from paper or foil or put the product in a sterilization container.

Respect requirements for sterilization packaging according to EN ISO 11607 and EN 868.

4.9. Sterilization
Only products that have been cleaned and disinfected are allowed to be sterilized.
Only apply steam sterilization in an autoclave (fractioned pre-vacuum with sufficient product drying) for this sterilization.

Adjust sterilization parameters:
• Minimum 134°C and maximum 137°C in saturated steam.
• Holding time at least 5 Min. until max. 20 Min.
• Drying in vacuum for at least 10 Min.
Example: Sterilizer Class B, manufacturer: Tuttnauer.
Respect data of the Sterilizer manufacturer regarding the load, handling, and drying times.

Do not apply hot air, EO-gas, Radiation, or Plasma for sterilization, or any other sterilization method for this product.
IMPORTANT:
Before usage, let the product cool to room temperature.

It is the sole responsibility of the user to maintain the sterile condition of the product after the sterilization process.
In case of the aforementioned chemicals and machines for cleaning, disinfection, or sterilization are not available, it is the responsibility of the user to validate the procedure actually applied.
Also, if a sterilization method other than that described above is applied, this deviating procedure has to be validated by the user accordingly.

4.10. Limitation of Reprocessing
Product life depends on wear and tear, handling, application time, damage and frequency and type of reprocessing.
Due to the materials and the intended use, this product may only undergo a maximum of up to 5 reprocessing cycles (see section "Visual and functional inspection").
Therefore, a visual and functional check must be carried out before each reuse.
Only an undamaged product may be reused.

5 Visual Inspection and Function Test
Before usage, the complete product, especially insulation and product tip have to be inspected for pressure points and damages.
Products exhibiting damage or pressure points are not allowed to be used and have to be replaced with new ones.
During and after application, tissue may adhere to the product, or sooting may be present on the distal end of the active electrode.
Such adhesions or sooting do not represent a reason for complaint and the product has to be exchanged for a new one.
Due to longer application time, mechanical forces, plasma seam, or similar, the product may exhibit deformation or abrasion of the insulation material.
Also, such aspects do not represent a reason for complaint and the product has to be replaced with a new one.
Blockage of the suction channel (if applicable), does not represent a reason for complaint.
A product, that exhibits a blocked suction channel (if applicable), has to be replaced with a new one.
Before usage, an electrical continuity test has to be done.
In case the product does not pass the electrical continuity test, it is not allowed to use the product any longer and has to be replaced with a new one.
If damaged, do not use these products.
Never lay this product on the patient or in their direct vicinity.

6 Exclusion of Repair and Modification
Defective products are not allowed to be repaired.
Unauthorized repair or repair is strictly prohibited and leads to an immediate loss of the manufacturer's warranty.
It is not allowed to bend products or manipulate the product in any other manner.
Especially products with a hook as an active part are never allowed to be bent.
This could lead to severe injuries to the patient, the user, or the third party.
The product is supplied non-sterile and must be cleaned, disinfected and sterilized before use according to the validated procedure (DIN EN ISO 17665) specified in these instructions for use.
Before usage, a visual inspection and function test must be carried out.
Only products that have passed the visual inspection and function test may be used again.
In particular, undamaged insulation and an undamaged electrode tip are decisive for whether a product may be used again or not.
See also label.

7 Packaging, Storage, Transport
The product must be stored in a clean and dry environment, they must be protected from direct sunlight.
Always handle the product with the utmost care when transporting, cleaning, disinfecting, maintaining, sterilizing, and storing it.
This is especially true for fine lines and other sensitive areas (e.g., insulation).
Do not store or transport the product together with sharp-edged or pointed objects.
Store only in protective containers with individual compartments or individual shrink-wrapped.
The maintenance of the sterile state after the sterilization process must be ensured by the operator.

8 Safety, installation and application
Only to be used by skilled medical professionals who have been introduced to the product.
In addition to the acknowledged benefits of HF technology, the application includes several risks that have to be attended to by the user.
Improper use and failure to follow the instructions for use may result in unintentional burns to the patient, user, or third parties.
Continued further education of the surgical personnel is recommended.
It is the operator's responsibility to ensure the training of users/surgeons in arthroscopic surgical techniques.
The HF generator, which is used in combination with the product, has to be checked for damages.
Follow the instructions outlined in the user manual of the manufacturer.

Preparation:
Prepare the HF generator according to the User Manual and place it for application. Proceed with the generator according to the User Manual.
Bring the sterilized product into the sterile field and confirm the integrity of the product's insulation by visual inspection.
Confirm the integrity of the insulation of the product by visual inspection.
At the beginning of the application, make sure that the cable of the handle is connected properly to the HF generator and that the correct settings are selected and displayed.

It is the responsibility of the user to select the lowest possible power setting on the HF generator to achieve the desired effect during the respective operation.
Note: when using products with suction, connect the suction adapter to a tube that leads to a suitable suction system (e.g., suction canister, wall suction system, etc.). The minimum vacuum target value is 300 mmHg.
This system is for the connect use of the HF neutral electrode, including patient protection and monitoring, monitoring of the HF neutral electrode and all other instructions, safety instructions and warnings in the instructions for use of the HF neutral electrode must be observed.
Make sure that there are no flammable substances (anesthetics, oxidizing gases, endogenous gases, etc.) in the immediate vicinity during the entire application, otherwise there is a risk of explosion.
Ensure correct patient positioning, i.e., using insulating OP drapes that are liquid-tight, dry, and absorbent.
Insulate conductive areas and contact points from the patient.
Place dry cellulose tissues in skin folds, breast folds, and between the extremities.
Any fluid that possibly pools in body cavities has to be removed before the procedure.
Only use non-flammable disinfectants.
Do not use alcohol-based lotions, for example.
Remove all body jewelry from the patient.
Putting a band-aid over the body jewelry is not sufficient!
All oxygen connections must be tight and leak-proof during the procedure.
As long as the product is not applied, place it on a dry, clean and non-conductive and well-visible surface, that is not in contact with the patient.
Never place the product on the patient.
Unintended activation of the product can lead to unintentional burns or other injuries to the patient, the user, or a third party.
Never wrap cables around the patient and never lay cables over the patient.

Operation:
Insert the product tip through the arthroscopic portal into the surgical field of vision. Ensure the tip of the product is fully submerged in a conductive fluid.
Only activate the product in the immediate vicinity of the tissue to be treated. **WARNING:** an active product tip must NEVER be outside of the user's field of vision during the surgery as this can cause unintentional tissue damage.
Whenever the product is not in use, place it on a dry, clean, non-conductive, and highly visible field that is not in contact with the patient.
Accidental activation of the product may cause burns to the patient.
Do not touch the tip of the instrument during the entire use.
No continuous activation.
Apply only short activation times.
At least up to 6 minutes cumulative ablation time and no longer than up to 6 seconds single activation.
Keep on longer breaks between activation phases.
Only adjust low power settings.
All electro-surgical instruments can potentially cause muscle stimulation during use.
The design of this product has been chosen to minimize the risk of this unwanted effect.
Nevertheless, muscle stimulation may cause unexpected patient movement in the surgical field.

End of surgery:
At the end of the surgery, remove the product carefully from the patient's body and confirm the completeness of the system.
Switch off the HF generator.
Remove the plug of the cable completely from the socket.

9 Manufacturer
REGER Medizintechnik GmbH
Gewerbestrasse 10 | 78667 Villingendorf | Germany.
Tel.: +49 (0) 741 270 698-0 | Fax: +49 (0) 741 270 698-10
info@regerr-med.de | www.regerr-med.de

10 Swiss Authorized Representative
CH REP
CHRN-AR-20001140
DR MEDICAL AG
Weissensteiner, Bau 6 | Ziegelmatteilweg 6 | 4503 Solothurn | Switzerland
Tel.: +41 (0)32 621 60 01 | Fax: +41 (0) 32 621 60 02
info@med-ag.ch | www.dr-medical.ch

11 Warranty
Our products satisfy the highest quality standards. Liability and warranty are excluded for all products that have been modified in any way, not applied according to their Intended Use, that has been handled or applied improperly, or in case of any other deviation from the instructions outlined in this IFU.
Furthermore, the manufacturer denies any liability for any accidental, intentional damage or damage or loss arising out of the use of the product. Additionally, all liability and warranty are excluded in case our product was repaired by a company that has not been authorized by us. Unauthorized repairs are strictly prohibited.

12 Return
Returned products will be accepted only if they are marked as "hygienically safe" or "not contaminated" and have been securely packaged for shipping.
Use our Return Form for returns.

13 Disposal
Disposal of the product, their packaging material as well as any accessory has to be done according to the applicable country specific requirements, regulations and laws.
Additionally, the applicable requirements of the respective clinical place(s) in regards to disposal of medical devices have to be followed as well.

14 Regulatory Remark
Due to regulatory reasons, we would like to inform users and patients, that serious incidents in connection with our medical device have to be reported to the manufacturer and competent authority of the EU Member State where the user and/or patient are residents.

15 About these Instructions for Use
Throughout the period of usage these IFUs must be kept freely accessible for the user.
For a current revision of these IFUs, please contact our customer service.
Changes reserved.
This medium is not intended for users in the USA.

©2022 REGER Medizintechnik GmbH All rights reserved. All company names, service marks, trademarks, product names and logos appearing on this document are the property of their respective owners and shall be treated accordingly.