
GEBRAUCHSANWEISUNG ARGON
INSTRUCTIONS FOR USE ARGON
MODE D'EMPLOI ARGON
INSTRUCCIONES DE USO ARGON
ISTRUZIONI PER L'USO ARGON
INSTRUKCJA KORZYSTANIA Z ARGONU

BOWA
E I N F A C H S I C H E R

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Flexible argon probes and connection cables for use with the HF surgical device and argon control device

Flexible argon probes		
TBS flexible argon probe (Ø 1.5 mm x l 1.5 m)	932-048	932-148
GIT flexible argon probe (Ø 2.3 mm x l 2.2 m)	932-049	932-149
GIT flexible argon probe (Ø 3.2 mm x l 2.2 m)	932-050	932-150
TBS flexible argon probe (Ø 2.3 mm x l 1.0 m)	932-118	932-154
GIT flexible argon probe (Ø 1.5 mm x l 3 m)	932-051	932-151
GIT flexible argon probe (Ø 2.3 mm x l 3 m)	932-052	932-152
GIT flexible argon probe (Ø 2.3 mm x l 2,2 m 90°)		932-153
Cable		
Connection cable for flexible argon probes	932-045	932-145 932-146

Liability

These devices comply with Directive 93/42 EEC on medical devices and are based only on original BOWA products.

Read the instructions for use carefully before using the devices.

Read all information on safety and precautionary measures before use!

Intended use

These devices may be used only in medical institutions by trained medical personnel.

The connection cable for flexible argon probes is a reusable, steam-sterilized supply line for flexible argon probes. The connection is made via an electrical line to the HF generator; the connection to the argon gas supply unit is made via an integrated hose adapter. The connection line with the flexible probes can be used in conjunction with HF generators for peak voltages up to a maximum of 5,200 Vp.

The device, together with the flexible probes, is designed for use in HF APC surgery. The argon control unit and the HF generator are activated using the foot switch.

DE

EN

FR

ES

IT

PL

Contra-indications

Do not cut metal stents.

Safety information

- Never apply the gas discharge jet of the flexible argon probe directly to the tissue during gas activation.
- ***There is a risk of gas emphysema!***
- The gas discharge jet should not be moved vertically over the tissue and activated. An angle of approximately 45°C works well.
- Do not blow argon into the vascular system (**risk of embolism!**)
- For endoscopic applications the probe tip must project out of the endoscope during activation.
- High HF voltages are required for gas ionization: you must ensure that the electrical insulation of the lines and applicators is intact.
- Attention! There is a risk of explosion and infection from endogenous gases. The patient must undergo intestinal cleaning before endoscopic procedures are carried out.
- Attention! Burn due to deflagration during bronchial use during ventilation with highly concentrated oxygen (100%). See study of HF surgery by BfArM 2009.
- Electrode wires that project distally can lead to mechanical tissue injuries and damage operating gloves.
- ***The argon plasma application is a monopolar HF surgical application; compliance with the safety rules for monopolar high-frequency surgery is required.***
- Argon of grade 4.8 must be used.
- The probe tips can be hot immediately after activation.
- ***There is a risk of burns!***
- Do not lay flexible probes on the patient.
- ***There is a risk of burns due to inadvertent HF activation!***
- Use only original BOWA accessories!

DE

EN

FR

ES

IT

PL

Instructions for Use

Flexible argon probes are used only with spray coagulation current subject to visual control!

Always operate the flexible argon probes using the blue pedal on the foot switch!

To „ignite“ the argon plasma, move the distal probe opening to within approximately 0.5 cm of the tissue. During the argon plasma coagulation procedure („beaming“), this distance can be lengthened to between 2 and 2.5 cm before the plasma is extinguished.

Use in flexible endoscopy

Use the flow rate given for the respective application as given in the table:

REF	Ø	FLOW l/min
932-048, 932-051, 932-148, 932-151	1,5 mm	0,2 - 0,5
932-049, 932-052, 932-118, 932-149, 932-152, 932-153, 932-154	2,3 mm	0,3 - 0,6
932-050, 932-150	3,2 mm	0,4 - 0,7

Connection of the flexible argon probes

Before beginning, ensure that the HF generator and the argon control unit are ready for operation. The gas tank valve must be opened and the gas supply sufficient.

1. Remove the connection cable and flexible argon probe from the container so that they remain sterile.
2. Check the instrument visually for damage to the insulation and the distal end of the flexible probes. Do not use damaged instruments! BOWA warns explicitly against repairing damaged sections. The instrument is operated at high coagulation voltages! Any brown discoloration at the distal end of the probe are due to the production process and does not constitute dirt or any other impairment of the product.
3. Connect the 3-pin plug with the HF surgical device. Generators with Plug'n Cut COMFORT can recognise BOWA COMFORT instruments and automatically select the appropriate parameters.
4. Insert the desired flexible argon probe with the plug into the coupling of

the connection cable. Ensure that the plug is pushed into the coupling until it hits the stop; a small gap (approximately 1 mm) will remain between the coupling and the plug.

5. Establish the argon gas connection to the argon control unit. To do this, connect the male Luer lock connector of the connection line to the female Luer lock gas discharge of the argon gas control unit. The use of a sterile filter (article no. 830-050) between the argon unit and the connection plug of the argon handle is recommended.
6. Activate the purge function of the argon gas control unit to flood the gas feed line of the connection cable with argon gas.
7. Activate the blue pedal of the foot switch and check to ensure that the system (activation of the argon gas control unit, generator and instruments) is working properly. Keep the tip of the probe at a distance of approximately 5 mm from any exposed, uninsulated metal part. It must not have any electrically conductive contact to the patient or other persons. Upon activation an argon plasma stream jet must form between the metal part and the tip.

The device is ready for operation.

Instrument preparation

The devices are not sterile! Clean and sterilize them before using them for the first time!

Non-sterile instruments must not be used without prior cleaning / disinfection and sterilization. Effective cleaning / disinfection is the crucial prerequisite for effective sterilization of the instruments.

As part of your responsibility for instrument sterility during use, you must ensure that only sufficiently validated processes for the specific devices and products are used in the cleaning/ disinfection and sterilization and that compliance with the validated parameters is ensured for each cycle. You must comply with the statutory regulations applicable in your country as well as the hygiene regulations of the hospital / clinic.

DE

EN

FR

ES

IT

PL

Overview:

1. Pre-cleaning

Soaking

The instruments must be soaked immediately after use (within 1 hour at the latest). Use only aldehyde-free disinfectants (to prevent fixation of contamination with blood) with tested effectiveness (e.g. with German Society for Hygiene and Microbiology (DGHM) or FDA approval and the CE mark) that are suitable for instrument disinfection and compatible with the instruments.

→ Recommended operating supplies for manual pre-cleaning, Page 24

The disinfectant used for soaking is only for the protection of persons and cannot replace the disinfection procedures to be carried out after pre-cleaning.

WARNING
Risk of infection

Manual pre-cleaning

Rinse off superficial impurities and coarse particles under running water, using plastic cleaning brushes if necessary.

Do not use scouring agents or metal brushes to clean the instruments.

Flush the inner adapter at least 3 times with approximately 50 ml (syringe) of cleaning fluid from the proximal to the distal end.

→ Recommended operating supplies for manual pre-cleaning, Page 24

DE

EN

FR

ES

IT

PL

After pre-cleaning, rinse the instruments thoroughly with fully desalinated or distilled water and flush the inner adapter 3-5 times with approximately 50 ml (syringe).

HF cable and Flexible argon probes not in the Ultrasound bath cleaning.

2. Cleaning / Disinfection

A machine process (disinfection unit) should be used for the cleaning and disinfection of the instruments. Manual procedures are not recommended due to their significantly lower effectiveness (use of ultrasound is not permitted!).

When selecting the disinfection unit, you must ensure,

- that the disinfection unit has tested effectiveness (e.g. with German Society for Hygiene and Microbiology (DGHM) or FDA approval and the CE mark in accordance with DIN EN ISO 15883),
- so that if possible a tested program is used for thermal disinfection (at least 5 minutes at 90 °C or A_0 value > 3000) (in the case of chemical disinfection, there is a risk of disinfectant residues on the instruments).
- that the program used is suitable for the instruments and includes sufficient rinsing cycles,
- that only sterile or low-bacterial-count (max. 10 bacteria/ml) and lowendotoxin (max. 0.25 endotoxin units/ml) water is used,
- that the air used for drying is filtered and
- that the disinfection unit is regularly serviced and inspected.

When selecting the cleaning agent system used, you must ensure,

- that this is fundamentally suitable for cleaning the instruments,
- that – insofar as no thermal disinfection is used – a suitable disinfectant with tested effectiveness (e.g. with German Society for Hygiene and Microbiology (DGHM) or FDA approval and the CE mark) is used and that this is compatible with the cleaning agents used, and
- that the chemicals used are compatible with the instruments

→ Recommended operating supplies for machine preparation, page 24

The concentrations and dwell times specified by the manufacturers of the cleaning and disinfecting agents must be complied with.


Procedure:

1. Place the instruments in the disinfection unit. Ensure that the cables are not kinked or pinched.

Connect the hollow cavities to flushing tubes and/or hoses to ensure sufficient rinsing of the instruments.


Use the **BOWA Luer rinsing adapter (REF 932-900 & 932-910)** as follows:

Step 1:

We recommend locking the probe into the adapter using slight pressure and then beginning the installation in the machine.

Step 2:

A secure connection of the adapter with the probe is facilitated in the washer via the lockable connection mechanism (Luer lock) of the two parts. The locking of the two parts is established by screwing the adapter into the washer connection.

DE

EN


FR

ES

IT

PL

Step 3:

2. Start the program.
3. Remove the instruments from the disinfection unit after the program has finished.
4. Inspect and pack the instruments immediately after removal (after additional post-drying in a clean place).

Because fluid penetrates the gas feed line of the connection cable as well as the probe during cleaning, blow out and dry the relevant adapters using filtered compressed air at a pressure of less than 3 bar after cleaning. To do this, place the compressed air attachment loosely on the Luer lock connector of the connection cable/probe.

Do not blow the air from the distal end to the proximal end!

3. Visual inspection and function testing

Check the devices..

→ Inspecting devices, page 17

4. Packaging

Prior to sterilization, the disassembled instruments must be packaged in suitable disposable sterilization packaging (single or double pouch) and/or a suitable sterilization container:

in accordance with DIN EN (ANSI AAMI) ISO 11607/
DIN EN 868-2...10 (previously DIN EN 868/ANSI AAMI
ISO 11607)

Suitable for steam sterilization (temperature resistance
up to 137 °C, sufficient steam permeability)
serviced regularly (sterilization container)

Sterilization in the shipping packaging (e.g. silicone protective cap, etc.) is prohibited!

DE

EN

FR

ES

IT

PL

5. Autoclaving

The instruments may be sterilized only after they are disassembled.

Only the sterilization methods listed below are to be used; other methods are not permitted.

Steam sterilization

- Fractionated vacuum method¹ (with sufficient device drying)
- Steam sterilizer in accordance with DIN EN 13060 or DIN EN 285
- Validated in accordance with DIN EN (ANSI AAMI) ISO 17665 (previously DIN EN 554/ANSI AAMI ISO 11134) (valid IQ/OQ (commissioning) and product-specific performance quality (PQ)).
- Maximum sterilization temperature 134°C (plus tolerance according to DIN EN (ANSI AAMI) ISO 17665 (previously DIN EN 554/ANSI AAMI ISO 11134))
- Sterilization time at least 20 minutes (at 121 °C) or 5 to 20 minutes at 132/134 °C.

The suitability of the devices for effective sterilization was proven by BOWA with the fractionated vacuum method using the sterilization times and temperatures above.

The basic suitability of the instruments for effective steam sterilization was documented by an independent accredited testing laboratory using the fractionated vacuum method (EuroSelectomat from MMM Münchener Medizin Mechanik GmbH in Planegg, Germany. The process described above was taken into account for this.

Hot air sterilization must not be used (because the instruments will be destroyed).

The manufacturer is not responsible for the use of other sterilization methods (such as ethylene oxide, formaldehyde, radiation and low temperature plasma sterilization). In this case the applicable standards are to be considered in each case (DIN EN ISO 14937/ANSI AAMI ISO 14937 and/or method-specific standards) and the suitability and basic effectiveness of the method (possibly including residue analyses of the sterilization medium) are to be documented taking into account the specific device geometry as part of the validation.

¹ The use of the less effective gravitation method must be backed up by an additional validation (longer sterilization times may be required).

6. Storage for use

7. Assembling parts/ Inspection / Function test in OR

Assemble the **disassembled** instruments as required for the relevant application. A visual inspection and function test must be carried out before each use.

→ Inspecting devices

→ Connecting the flexible probes, Page 17 and 18

Material compatibility / reusability

With appropriate care and inspection, the flexible argon probes and the connection cable can be re-used up to 20 times, the cable 200 times, the COMFORT cable 100 times (in each case with reference to a sterilization time of 20 minutes and a sterilization temperature of 134 °C); every use beyond that is at the user's own responsibility.

The service life of the devices was validated with the use of the recommended cleaning and preparation procedures. Improper handling and unsuitable preparation media reduce the service life.

Defective probes and cables may not be repaired. Always replace a defective part with a new one.

The manufacturer is not liable if these instructions are disregarded.

Recommended operating supplies

Cleaning agents and disinfectants for manual pre-cleaning:

Manufacturer	Product
Schülke & Mayr Robert-Koch-Str. 2 22851 Norderstedt, Germany	Gigasept Instru AF

The suitability of the devices for effective ultrasonic cleaning (5 min.) was documented by BOWA using an aldehyde-free combination cleaning and disinfecting agent (**Gigasept® Instru AF**).

DE

EN

FR

ES

IT

PL

Cleaning and disinfecting agent for machine preparation:

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Product

neodisher MediClean forte

BOWA recommends the use of neutral to slightly alkaline cleaning or cleaning/disinfecting agents without hazardous ingredients. Depending on the concentration, alcohol and/or aldehyde ingredient materials are permitted.

The suitability of the devices for effective cleaning/disinfection was documented by BOWA for the machine process (90°C, 5 minutes) using an alkaline cleaning agent with added tenside (**neodisher® MediClean forte**)

This applies in the following conditions: The described cleaning and disinfecting method is applied when cleaning/disinfecting. Cleaning agents must offer the proven effectiveness and be compatible with the system (disinfectant and cleaning agent) used in the specified test.

When alkaline cleaning agents (pH 9.5 – 11.5) are used, discoloration of metallic surfaces may occur. This does not affect the proper functioning of the device. Do not use strongly alkaline cleaning agents (pH > 11.5).

Storage and disposal

Storage:

Store dry at room temperature.

Do not expose the equipment to strong mechanical effects such as impacts, dropping or striking.

Storage conditions:

Store dry at room temperature.

(max. temperature -25°C to +50°C, relative humidity 0 to 75% and air pressure of 500 to 1600 mbar)

Keep out of direct sunlight.

Avoid continuous exposure to X-rays.

DE

EN


FR

ES

IT

PL

Duration of storage:

The storage period for sterilized instruments depends on the type of packaging and the storage conditions.

Medical devices, packaging material and accessories must be disposed of according to the applicable domestic laws and regulations.

DE

EN

BOWA-electronic GmbH & Co. KG
Heinrich-Hertz Straße 4-10
D-72810 Gomaringen

FR

You can obtain more information under the following number:
+ 49 (0) 70 72-60 02-0

ES

IT

or on the Internet at: www.bowa-medical.com

PL

Marketed through the authorised specialist medical dealer

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Bowa-electronic GmbH & Co. KG