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
The manufacturer and vendor of this fiberoptic light cable does not accept any liability for direct or consequential damages caused by the improper use or handling, in particular as a result of the failure to follow the instructions below or by improper processing and maintenance.

	The high light concentration at the end of the fiberoptic light cable generates heat. Never place the end of the fiberoptic light cable onto the cover or the skin of the patient, because the light intensity may cause burns to the patient's skin or set the cover on fire.
	Fiberoptic light cables are supplied in non-sterile condition. Therefore, thorough cleaning, disinfection and/or sterilisation are required before the first use as well as before any subsequent utilisation.
	Fiberoptic light cables must not be cleaned in the ultrasonic bath.
	Damage of the coating causes damage to the fiberoptic light cable! Strong kinking causes glass-fiber-breakage and destroys the fiberoptic light cable.
	Damages of the fiberoptic light cables arising as a result of improper handling are not covered by the warranty. It is essential to avoid hitting in particular the adapters and bending the light cables, because this may damage the coating and the optical fibres.
	The cross-sections of the cables must match the diameter of the light input of the endoscope to prevent damages at the light output of the cable caused by reflection. 3,5mm - thin endoscopes, such as e.g. for arthroscopies, urological imaging, sinuscopies, etc. 4,9mm - thick endoscopes, such as e.g. for laparoscopies, etc. 

Handling and Application / Indication

Fiberoptic light cables are used to transmit light energy from a cold light source to an endoscope. (Area of application exclusively outside the body). The medical device may only be used by physicians and medical assistant personnel who are appropriately qualified. Indications: The products have no medical indication. Contraindications: There are no other known contraindications directly related to the medical device if the instructions for use are observed.

Guideline conformance

This medical device bears the CE mark pursuant to the MDR 2017/745. Furthermore, every fiberoptic light cable has a separate serial number. Any adaptations supplied are variable components of cold light cables and do not carry a CE mark.

Combination

The RfQ Light Cables are designed for use with halogen or xenon / LED light sources. It is important to note that glued fiber optic cables may only be used in halogen light sources. The following combinations are used.

RfQ-Fiber Light Cable	Connection Light Input	Connection Light Output
REF 00-2xxxx-xx	ACM/Aesculap/Dyonics/ Heine/HSW/Medicon/Olympus/ Pilling/RFQ/Schölly/Schott/ Storz/Stryker/Volpi/Winter&IBE/ Wolf/Zeiss	ACM/Aesculap/Dyonics/ Heine/HSW/Olympus/RFQ Storz/Stryker/Wolf/Zeiss

Function test

Procedure:

- Hold one side of the fiberoptic light cable pointing toward a light source (e.g. window). For fiberoptic light cables, check the number of dark dots at the other end. These dark dots indicate broken optical fibers of the fibre optic bundle. Several broken optical fibres will not result in a noticeable impairment of the image quality. However, starting from about 20 to 25% of the light output surface, the light transmission and hence the image quality start to deteriorate considerably.
- Check the fiberoptic light cable and the cable ends for damages such as burns or erosions.

Information:

- Please note that the unscrewable adapters may potentially be re-used when ordering a new fiberoptic light cable.


Cleaning preparation

To achieve an optimal cleaning and disinfection outcome, both adapters should be disassembled from the fiberoptic light cable

Procedure Disassembly:

- To disassemble, screw the adapters from the respective end of the cable by turning anti clockwise.

Procedure Assembly:

- To assemble, screw the adapters onto the respective end of the cable by turning clockwise.

Information:

- Place the fiberoptic light cables onto a disposal screen in such a way that no rinsing-related shadows are created!
- Remove major contaminations from the fiberoptic light cable at the operating table.

Manufacturer Specifications for the Preparation Procedure of Medical Devices According to DIN EN ISO 17664:2021 (Preparation Instructions)

These recommended and validated preparation instructions primarily apply to the reusable universal cold-light cable of endoscopy, which is referred to as 'product' in the following.

	The use of Xenon- or LED- cold light sources is not permitted for the RfQ light guide cable (= LLK) bonded on the light entry, and furthermore, this type of light guide cable must not be processed with thermolabile sterilization processes (plasma).
	The following procedure may not be used for the preparation of the products: Ultrasonic cleaning, flash autoclaving and hot air sterilisation.
	Limitation of re-preparation: Frequent preparation has little impact on these products. The end of the product's life is normally caused by wear and damage from use. (-> see also further -> point reusability)

Recommended & validated preparation instructions for the product¹

Preparation at the location of use

Wipe off surface contamination with a non-shedding single-use towel. Properly put away the product in order to prevent damages.

Storage and transport

It is recommended to store the product in containers which are suitable for the transport. The transport into the preparation rooms can take place wet as well as dry. For a dry disposal, make sure that no residuals dry on the surface. Close the container. The cleaning should take place within three hours. For a wet disposal, the cleaning must take place within one hour, and the recommended combined cleaner and disinfectant (see manual cleaning) must be used. Do not exceed the time windows stated above. The product must not be kinked or bent in a radius that is too narrow. It is very important to prevent bumping, especially of the fused light entry, as this can lead to the destruction of the fusing.

Cleaning preparation

Remove the adapters from the product, as a separate, manual cleaning and disinfection of the individual parts will take place. Properly put away the individual parts in order to prevent damages.

Precleaning

Required aids: Sterilization tray, water bath, Cleaning solution: Neodosher Medizym 0,5% (v/v), Cold Tap water (15-20°C, max. 45°C), Soft nylon brush.
Procedure: **1.** Place the individual parts into the sterilization tray and then transfer them into an immersion bath with the cleaning solution of Neodosher Medizym 0,5% (v/v). **2.** After 5 minutes soaking time: **a.** Brush off under the water surface using a soft nylon brush. **b.** Dry the exterior with a non-shedding single use towel or swab. **c.** Rinse with cold tap water.
Do not use any metal brushes or cotton wool applicators made from metal for the manual cleaning. Do not use any other instruments for cleaning the optical surfaces. After the manual Precleaning, inspect all the individual parts for damages (see checking and maintenance).

Machine cleaning and Disinfection

The international standards series EN ISO 15883 and national guidelines state that only validated machine-based cleaning and disinfecting procedures should be used. The requirements imposed upon washers and disinfectors are described in the EN ISO 15883 series of standards.

Required aids:

- Washer-disinfector unit e.g. Miele PG 8535
- Cleaning solution: Neodosher Mediclean 0,5 % (v/v)
- Cold Tap water
- Cold DI water

Procedure:

1. Pre-rinsing 1 min. with Cold Tap water
2. Soaking time 1 min.
3. Cleaning with Neodosher Mediclean 0,5 % (v/v) for 600 sec. at 55°C (131°F)
4. Post rinsing 1 min. with Cold DI water.
5. Soaking time 2 min.
6. Thermal disinfection at 90°C (194°F) for 1 min. (= A0 600)
7. Drying 20 min. at 100°C (212°F)

¹ According to the test reports 5561 and -5563 from 28/08/2018 from the CleanControlling Medical GmbH & Co. KG – 78576 Liptingen / Germany

Alternative the 'Vario TD programme' for the machine cleaning and disinfection procedure may be used.

Excerpt from <http://www.miele-professional.de>

„After an intensive pre-rinse with cold water, the cleaning phase takes place at a temperature of up to 55°C – five minutes hold time. The process concludes with the thermal disinfection with 90°C and at least 60 seconds (1 minute) temperature hold time. For optimal care of the instruments, the final rinse should preferably take place with DI water without rinsing agents.“

Note

The universal light guide cable should be stored in a suitable container (sterilization tray / basket) in the machine, in order prevent damages of the product. Corrosion on the product can occur from using chloride-containing water; this is why the final rinse procedure should take place with fully demineralised water. The operating instructions and loading instructions of the manufacturer must be precisely adhered to during the cleaning and the following thermal disinfection. The cleaning agents used must be exactly dispensed according to the prescription of the manufacturer. The disinfection temperature must not exceed 93°C.

Drying

Required aids: Non-shedding towel, Compressed air

Procedure: **1.** Put the cleaned and disinfected individual parts back together. **2.** Dry the product by means of the non-shedding towel and/or the compressed air. **3.** Properly put away the parts in order to prevent damages.

Checking & Maintenance

We recommend a visual inspection in which the product should be checked for the following defects: Damage, sharp edges, loose or missing parts, rough surfaces, residues of cleaning agents and disinfectants (residues must be removed), labels and markings that are safe and proper Use required must be legible.

Beware of damaged and incomplete product (injuries of the patient, user or third parties are possible). Perform checks before and after each use. Do not use product that is damaged and / or incomplete or contains loose parts. Send damaged product with loose parts for repair. Do not make your own repair attempts.

We also recommend a visual inspection of the optical surfaces on:

1. light output, **2.** broken fibers (appear as black dots on the cold light connection) -> the proportion should not exceed 30%, **3.** coverings on the glass surfaces (may worsen the light transmission) -> clean glass surfaces (see manual cleaning & disinfection)

Note Do not carry out any repair attempts yourself if coverings can not be removed with the recommended cleaning agents and disinfectants (see Cleaning & Disinfection). Regular cleaning with the 70% alcohol solution after each preparation prevents deposits

We recommend the following maintenance instructions before sterilization: Reassemble cleaned and disinfected parts.

Manual cleaning. Not recommended!

Required aids: Sterilization tray, water bath, Cleaning solution with disinfecting effect: Sekusept 4%, Cold Tap water (15-20°C, max. 45°C), DI (fully de-mineralized) water, Non-shedding single-use towel or swab

Procedure: **1.** Thoroughly rinse off the individual parts with tap water (max.45°C). **2.** Place the individual parts into the sterilization tray and then transfer them into an immersion bath with the self-acting disinfecting cleaning solution of Sekusept 4%. **3.** After the recommended soaking time according to the manufacturer specifications of the cleaning solution: **a.** Rinse each universal cold-light cable for five minutes with DI water. **b.** Dry the exterior with a non-shedding single use towel or swab. **c.** Clean the mechanical parts and optical surfaces (= light entry and light exit) with a soft towel or soaked cotton ball and the 70% alcohol solution, provided that a disinfection does not take place after the cleaning. **4.** Properly put away the individual parts in order to prevent damages.

Note Do not use any metal brushes or cotton wool applicators made from metal for the manual cleaning. Do not use any other instruments for cleaning the optical surfaces. After the manual cleaning, inspect all the individual parts for damages (see checking and maintenance).

Manual disinfection. Not recommended!

Required aids: Sterilization tray, disinfection bath, Disinfection solution, 70% alcohol solution (ethanol, isopropyl), DI (fully de-mineralized) water, Non-shedding single-use towel or swab.

Procedure: **1.** Place the individual parts into the sterilization tray and then transfer them into an immersion bath with disinfection solution. (Refer to the specifications of the chemical manufacturer for the concentration and soaking time of the deployed disinfectant.)

2. Then thoroughly rinse the universal cold-light cable for five minutes with DI water.

3. Dry the exterior with a non-shedding single use towel or swab. **4.** Clean the mechanical parts and optical surfaces (= light entry and light exit) with a soft towel or soaked cotton ball and the 70% alcohol solution. **5.** Properly put away the individual parts in order to prevent damages.

Note Disinfectants which contain acetic acid or chlorine compounds may not be used. After the manual disinfection, inspect all the individual parts for damages (see checking and maintenance). Observe the manufacturer's specifications for the disinfectant regarding the: - Disinfectant effectiveness, - Concentration, - Soaking time and -Product life

Sterilisation

Required aids: Storage system for the sterilisation, Steam sterilizer in accordance with DIN EN 285, Sterile packaging

Procedure: The medical device must be sterilised in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) at 132°C – max 138°C with a minimum exposure time of 2 minutes and a maximum of 5 minutes.

Note The cleaned and disinfected individual parts must be put back together before the sterilisation. The product must be sufficiently clean and dry. Observe the manufacturer's specifications / instructions for use of the sterilizer (maintenance intervals must be adhered to). Observe the national and international standards for the sterile packaging. Brand-new light guide cables must be sterilised before their first use.

Storage

After disinfection, store the product according to the following conditions:

Completely dry, Protected from dust, In a closed container, In a germ-free environment

Note If the product has been stored for several days, it must be disinfected again before the sterilisation.

After sterilisation, store the product in sterile goods packaging as follows:

- Protected from humidity and temperature fluctuations, -Protected from direct sunlight,
- Protected from dust

Note Improper storage can lead to loss of sterilisation - the manufacturer assumes no liability in this case.

Reusability

The service life of RfQ cold-light cables is largely determined by wear or damage from use. The number of reprocessing cycles is not quantitatively limited. Before reuse, the specified functional test must be carried out. Reprocessing in a validated reprocessing procedure is assumed here. In order to generate an orientation value for a quantifiable service life, the reference cold-light cable was reprocessed 100 times by the manufacturer according to the specified reprocessing procedure without loss of functionality.

Additional information

In order to prevent damages while transporting the products, we recommend using the original packaging for sending. The above listed instructions were validated by a medical device manufacturer for the preparation of a medical device whose reuse is validated as SUITABLE. It is the responsibility of the preparer to ensure that the equipment, materials, and personnel in the preparation facility, which are used with the actual carried out preparation, achieve the desired results.

This requires validation and routine monitoring of the procedure. Additionally, each deviation from the provided instructions should be carefully evaluated by the preparer for effectiveness and possible negative consequences.

Annex: Further possible preparation procedures

Reusable universal light guide cable					
Procedure			Universal light guide cable	Adapter	
Preparation at the location of use	State		Wet disposal	●	●
			Dry disposal	●	●
Decontamination	Preparation		Selection of the cleaning and disinfection procedure and removal		
	Cleaning	Procedure	Manual not recommended	●	●
			Machine 01)	●	●
			Ultrasonic	○	○
	Cleaner		Alkaline not recommended	●	●
			Citric acid not recommended	●	●
			pH neutral / enzymatic	●	●
			Containing acetic acid	○	○
	Rinse	Procedure	Manual not recommended	●	●
			Machine 01)	●	●
	Disinfect 03)	Procedure	Manual not recommended	●	●
			Machine 01) - Thermally max. 93°C	●	●
	Dry	Procedure	Manual not recommended	●	●
			Machine 01) - Tmax in 100°C	●	●
Sterilisation	Steam	Fractionated vacuum process 132°C (max.) 138°C / 273°F 01)		●	●
	Low-temperatur 04) not recommended	Formaldehyde (FA)		●	●
		Ethylene oxide (EO)		●	●
		STERIS®		○	○
		STERRAD®		○	○
01) See the recommended product preparation, pages 1-4			Legend ● permitted ○ not permitted		
02) DI (fully de-mineralised water) recommended for the final rinse					
03) Released disinfectant component groups: Aldehyde, QAC's, Alcohol, Amine derivatives					
04) under own responsibility (without present validation), however the reference: no recommendation of mutual employment of STERIS® and STERRAD®					

Disclaimer

The manufacturer of the fiberoptic light cable does not accept any liability for direct or consequential damages caused by the improper use, handling or improper processing, sterilisation and maintenance. Failure to follow the instructions, improper handling or use of products supplied by the manufacturer of the fiberoptic light cable for unintended purposes will completely void the warranty. The manufacturer of the fiberoptic light cable does not accept any liability whatsoever for damages arising hereof. All serious incidents must be reported to the manufacturer (vorkommnis@rfq.de) and the competent authority. A "serious incident" includes incidents that directly or indirectly had, could have had, or may have had any of the following consequences: A serious public health hazard, Temporary or permanent serious deterioration in the health of a patient, user or third party, Death of a patient, user or third party. Products must be decontaminated before disposal.

Status of Information

These operating instructions were issued in accordance with up-to-date facts and the date of revision is: **02/2023**

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RfQ-Medizintechnik GmbH & Co. KG