

Our ^{new} mediPlogs



New Design – Improved Quality!

Our feet carry us 130,000 kilometers through life, almost three times around the globe. They lift 2,500 tons per day – for an average of 78 years. With 26 bones, 30 joints, 60 muscles, 100 ligaments and 200 tendons, they are a true masterpiece of the human body. Are you impressed? Us too! medimex stands for decades of experience in developing surgical shoes of the highest quality.



Their good sole profile ensures safe posture.



The special ventilation system ensures that your feet stay fresh.



The adjustable heel strap keeps your foot securely in the shoe.



The orthopaedic insole supports freedom from fatigue and provides relief.

mediPlogs: Our classic shoes - always impressive-

A special ventilation system keeps your feet cool

- ensures a pleasant environment for your feet
- the characteristics of the lateral ventilation holes prevent extensive penetration of liquids into the shoe

Heel lift and orthopaedic insole

- ensure fatigue-free walking and standing
- relieve the musculoskeletal system at the operating table

Machine washable up to 60 °C

- theatre shoes made of PU (polyurethane) are usually washed chemothermally at 60 °C with 1% disinfectant solution and a soaking time of 5 minutes

Made of anti-static PU

- foamed without CFC
- provides safety and protects the environment

Complies with applicable essential health and safety requirements

- Corresponds to the requirements of the harmonized standard EN ISO 20347:2012
- PPE protection class: OB SRA A E
- PPE risk category: II

Labeling service for mediPlogs

- on request the **mediPlogs** can be labeled with name/department abbreviation letters

Colour	Size	Pair/ packing unit	Item no. without heel strap*	Item no. with heel strap*
green 	35–48	1	8M0374__	8M0384__
blue 	35–48	1	8M0375__	8M0385__
white 	35–48	1	8M0378__	8M0388__
fuchsia 	35–48	1	8M0376__	8M0386__
light blue 	35–48	1	8M0377__	8M0387__

*Please include the size with the item number when ordering.



medimex 

medimex GmbH

In den Fritzenstücker 9–11, D-65549 Limburg
fon +49.6431.7302 300, fax +49.6431.7302 235
www.medimex.de