



Kabelmatrix Anschlusskabel

Geräteseitiger Anschluss			BOWA	Conmed	Erbe VIO	Erbe ICC	Erbe ACC	Erbe T-Serie	Martin	Valleylab
 Monopolarer Handgriff	Standard	 = Ø5,0 mm	•*	—	•	•	•	—	—	—
	3-Pin-Anschluss		•	—	(int.)	(int.)	(int.)	—	—	—
	M-Anschluss	 12,5 mm 4,0 mm	—	—	—	—	—	•	—	—
	T-Anschluss	 = Ø4,0 mm	—	—	—	—	•	—	—	—
 Bipolare Pinzetten	Standard	 12,5 mm 4,0 mm	•*	—	•	•**	•**	•	—	—
	2-Pin-Anschluss (22 mm)	 22,0 mm	—	—	(int.)	—	—	—	—	—
	2-Pin-Anschluss (28 mm)	 28,0 mm	•	•	—	—	—	—	(int.)	•
	M-Anschluss	 6,0 mm 2,0 mm	—	—	—	—	—	•	—	—
 Bipolare Arthroscopie-Elektroden	Standard	 = Ø6,3 mm	—	—	•	•	•	•	•	—
	International	 10,0 mm 2,5 mm	•	•	(int.)	(int.)	(int.)	—	—	•
	Standard	 12,5 mm 4,0 mm	•*	—	•	•	•	•	—	—
	2-Pin-Anschluss (28 mm)	 28,0 mm	•	•	—	—	—	—	(int.)	•
 Neutral-elektroden	M-Anschluss	 6,0 mm 2,0 mm	—	—	—	—	—	•	—	—
 Bipolare Scheren	Standard	 = 5,0 mm	•*	—	•	•	•	—	—	—
	Bovie	 = 8,0 mm	•	•	(int.)	(int.)	(int.)	(int.)	(int.)	•
	4 mm Buchse	 = 4,0 mm	—	—	—	—	—	•	—	—
	Standard	 = Ø 5,0 mm	•*	—	•	•	•	—	—	—
 MLC-Instrumente	Standard	 = Ø 5,0 mm	•*	—	•	•	•	—	—	—
	Standard	 = Ø 5,0 mm	•*	—	•	•	•	—	—	—
	Standard	 = Ø 5,0 mm	•*	—	•	•	•	—	—	—
	Standard	 = Ø 5,0 mm	•*	—	•	•	•	—	—	—
 MTW-Schlingen	Standard	 = Ø 5,0 mm	•*	—	•	•	•	—	—	—
	Standard	 = Ø 5,0 mm	•*	—	•	•	•	—	—	—
	Standard	 = Ø 5,0 mm	•*	—	•	•	•	—	—	—
	Standard	 = Ø 5,0 mm	•*	—	•	•	•	—	—	—
 Olympus-Schlingen	Standard	 = Ø 5,0 mm	•*	—	•	•	•	—	—	—
	Standard	 = Ø 5,0 mm	•*	—	•	•	•	—	—	—
	Standard	 = Ø 5,0 mm	•*	—	•	•	•	—	—	—
	Standard	 = Ø 5,0 mm	•*	—	•	•	•	—	—	—

- passender Anschluss
- * ARC 400 + 350
- ** 8M960057 ist nicht kompatibel
- (int.) international
- nicht kompatibel



medimex 

medimex GmbH

In den Fritzenstücker 9–11, D-65549 Limburg
fon +49.6431.7302 300, fax +49.6431.7302 235
www.medimex.de